



## NEW PRODUCTS



### AN-PD2 LCD Annunciator Display **NEW!**

Austco's new high resolution, LCD Annunciator provides audio and visual notification of active alarms with accompanying alert tones based on priority. The display supports multiple languages including the display of non-western characters. The AN-PD2 display allows alert messages, colors and tones to be completely customised and configured to suit your facility's needs.

The display is available as a single-sided unit or double-sided unit ideal for use in hallways, corridors and any space that requires viewing from both sides. If no active alarms are present, the display will default to a screen saver which can be an image, customizable message, date, or clock with customer branding.

The AN-PD2 plays clear alarm tones and other audio at up to 70 decibels at 1 meter with easily adjustable volume to assist with overall facility ambience. The non-reflective screen of the AN-PD2 helps eliminate glare and wide-viewing angle makes the display easy to read from a distance.

[DOWNLOAD CUT SHEET](#)




### IP-TSNS-2 Master Station **NEW!**

The new IP-TSNS-2 touch screen master station has a modern design with exceptional performance. It is designed for use with Austco's Tacera nurse call systems and is ideal for use in hospitals, aged-care facilities, treatment and outpatient centers, or anywhere requiring alarm management and communications.

The Android-based device runs Austco's Pulse Care software for intuitive alarm management, workflow, and communications.

The IP-TSNS-2 has a small footprint on a desktop, can be wall mounted and plug and play deployment for an easy out-of-the-box experience.

[DOWNLOAD CUT SHEET](#)



Contact your local Austco representative for more information on new products.



## Austco's Built-in RTLS


Austco's Built-in RTLS provides real-time locating, built into the Tacera nurse call system without the need for 3rd party hardware. Your nurse call hardware is also your RTLS infrastructure. With each interaction of badge and sensor, data is sent to the RTLS server, and Tacera records the badge location and time.



## Streamline clinical communication processes with workflows

Clinical workflow is one of the most powerful additions to a Nurse Call System. Austco's workflow features allow you complete control over your processes by allowing you to:

- Accurately monitor rounding compliance
- Streamline your communication channels
- Record pain checks and other care initiatives
- Push data to a patient's chart with a single button
- Indicate statuses like staff presence or patient falls risk by turning on the over-door light



## Have an idea for a News Story?

If you have a news story to share, we want to know about it! Please email your stories to [marketing@austco.com](mailto:marketing@austco.com) along with any photos, videos or other content related to your story.

## How to improve infection control through Austco technology

As infectious viruses such as the Coronavirus (COVID-19) continue to spread and impact healthcare facilities around the world, it is crucial to choose technologies that can help your organization combat the spread of infection.

Austco has designed its Nurse Call solutions to help with infection control in healthcare facilities. Austco Nurse Call Systems have:

- Antibacterial silicon buttons
- Spray and wipe maintenance
- Dip-sterilizable feature
- Easy maintenance and cleaning



Patient stations incorporate a cleaning mode for easy maintenance. Cleaning time and date is logged providing an important audit trail for infection control compliance.

Austco Dip-Sterilizable pendants have an IP67 rating ensuring the pendants are waterproof. The soft touch silicone case and buttons also have an anti-bacterial additive.

Vinyl washroom pull cords have a wipe-clean surface to aid in improved infection control.

Call point buttons are manufactured with an anti-bacterial additive with a spray-and-wipe clean surface.

## Austco's Sales Team is growing in the USA

We're pleased to welcome the newest members of our US Sales Team:



**Adam Kuehn**  
Direct Sales Rep

Adam is a seasoned sales professional, hailing from the telecommunications world with a zest for finding the right messaging to meet prospects where they are and a passion for crossing the finish line first!



**Annie Tamburro**  
Direct Sales Rep

Over the past 15 years, Annie has built and solidified her reputation as an accomplished sales and business development executive. Her sales experience covers industries from software, government, and artificial intelligence (AI).



**Tate Baird**  
Direct Sales Rep

Tate is a sales and business leader with 15+ years of Software, Telecom, IT, Marketing, Sales, and Business Development experience. His background also includes 7 years selling SaaS to mid to enterprise businesses.